

Biosafety Components

SNL Biosecurity Team
National Workshop on Biosecurity
Kuala Lumpur, Malaysia
1 June 2005







Biosafety Level 1

- Suitable for work involving well-characterized agents not known to cause disease in healthy adult humans and of minimal potential hazard to laboratory personnel and the environment.
- Examples:
 - Bacillus subtilis
 - Naegleria gruberi
 - Infectious canine hepatitis virus
 - E. coli





Biosafety Level 1: Facility Design

Requirements:

- Laboratories have doors
- Sink for hand washing
- Work surfaces easily cleaned
- Bench tops are impervious to water
- Sturdy furniture
- Windows fitted with flyscreens


Easily cleaned and decontaminated







Biosafety Level 1: Facility Design







Biosafety Level 2

- Suitable for work involving agents of moderate potential hazard to personnel and the environment
- Examples:
 - Measles virus
 - Salmonellae
 - Toxoplasma species
 - Hepatitis B virus





Biosafety Level 2: Facility Design

Requirements:

- Laboratories have lockable doors
- Sink for hand washing
- Work surfaces easily cleaned
- Bench tops are impervious to water
- Sturdy furniture
- Biological safety cabinets installed as needed
- Adequate illumination
- Eyewash readily available
- Windows fitted with flyscreens
- Location separated from public areas
- Ventilation directional
 - Air flows into lab without re-circulation to non-lab areas




Restricted access when work in progress





Biosafety Level 2: Facility Design







Biosafety Level 3

- Suitable for work with infectious agents which may cause serious or potentially lethal disease as a result of exposure by the inhalation route.
- Examples:
 - Mycobacterium tuberculosis
 - St. Louis encephalitis virus
 - Coxiella burnetii






Biosafety Level 3: Facility Design

Requirements:

- BSL-1 and 2 Facilities PLUS
 - Enclosures for aerosol generating equipment
 - Room penetrations sealed
 - Walls, floors and ceilings are water resistant for easy cleaning
 - BSC class II or III to manipulate infectious material






Respiratory protection may be appropriate





Biosafety Level 3: Facility Design







Summary of Biosafety Level Requirements

Table 3. Summary of biosafety level requirements

	BIOSAFETY LEVEL			
	1	2	3	4
Isolation ^a of laboratory	No	No	Yes	Yes
Room sealable for decontamination Ventilation:	No	No	Yes	Yes
 inward airflow 	No	Desirable	Yes	Yes
 controlled ventilating system 	No	Desirable	Yes	Yes
 HEPA-filtered air exhaust 	No	No	Yes/No ^b	Yes
Double-door entry	No	No	Yes	Yes
Airlock	No	No	No	Yes
Airlock with shower	No	No	No	Yes
Anteroom	No	No	Yes	_
Anteroom with shower	No	No	Yes/No ^c	No
Effluent treatment Autoclave:	No	No	Yes/No ^c	Yes
— on site	No	Desirable	Yes	Yes
 in laboratory room 	No	No	Desirable	Yes
 double-ended 	No	No	Desirable	Yes
Biological safety cabinets	No	Desirable	Yes	Yes
Personnel safety monitoring capability	No	No	Desirable	Yes

a Environmental and functional isolation from general traffic.

From: WHO LBM 3rd edition



b Dependent on location of exhaust (see Chapter 4).

Dependent on agent(s) used in the laboratory.

d For example, window, closed-circuit television, two-way communication.